

MEASUREMENT
FORM



Body Suit

Body Suit

Order/Reorder Form

All measurements taken with no tension and in centimetres only.


Design options



Female bra cup
with full scoop



Male pouch



Snaplock opening
- body brief



Postural support
- leg abduction




Uplift and double
abdominal panel



Open crotch



Velcro[®] opening
- body brief



Motifs

Design choice definitions

Full neck scoop: This option will increase the opening at the front, back and sides of the neck as well as decreasing compression over these areas.

Shoulder bracing: A diagonal pull to help extend the shoulders backwards and downwards to assist in improving posture as well as garment fit over the anterior chest wall.

Uplift panel: A 'v' shaped double layer of powernet extending

from the middle of the abdomen to assist in supporting a pendulous abdomen.

Double abdominal panel: A double layer of powernet over the abdomen to provide increased abdomen support.

Soft abdominal panel: A section of satin lining across the abdomen for pregnancy.

Measuring tips

Neck circumference (1): Measured midway between the head and the shoulders.

Scoop neck back: Measure from C7 to depth of scoop.

Scoop neck front: Measure from sternal notch to depth of scoop.

NB: Scoop neck. This option is not recommended when treating the upper torso as a large scoop will significantly reduce the pressure around the shoulder, upper chest and back.

V neck front: Measure depth required from sternal notch to point of V.

Attached powernet turtleneck: Measure height from sternal notch.

Shoulder circumference (2 & 3): With patient standing with hand on hip, place tape measure on top of shoulder directly above axilla, and pass through under the arm. Then lower arm down and continue tape to shoulder.

Arm circumference: Using Moss Jobbskin arm tape.

Chest circumference (4): Measured at the level of the axilla with arms placed down to the side.

Nipple line (5): For women please provide bra cup size.

Circumference over nipple line should be measured with patient wearing a bra.

Side seam to side seam (19 & 20): A horizontal contoured measure originating directly below the midpoint of one axilla across the back to the midpoint under the other axilla. (19) is positioned of chest circumference (4), whilst (20) is positioned at the height of the waist.

Leg circumferences: Begin with the patient lying down, knee straight and leg supported in elevation. Proceed using an Teal Jobskin leg tape. For an above knee garment, start tape at the gluteal fold. pleating only at the distal end if required. For a below knee garment secure your first tape at the point you wish the garment to end (generally ankle 1.5" tape, or base of calf 4.5" tape); pleating at the gluteal fold if required. Please indicate which tape falls closest to the popliteal crease.

Email, fax or return forms to:

**Smith & Nephew Pty Ltd
Australia**
Unit 1, Level 1, 100 Station Street
Nunawading
Victoria 3131, Australia
T 1800 772 577
F 1800 773 577
www.smith-nephew.com/australia
E Jobskin.cs@smith-nephew.com

**Smith & Nephew Ltd
New Zealand**
621 Rosebank Road, Avondale
PO Box 422, Shortland Street
Auckland 1140, New Zealand
T 0800 765 627
F 0800 465 278
www.smith-nephew.com/new-zealand

Garment returns

Smith & Nephew Pty Ltd Australia
Unit 1, Level 1, 100 Station Street
Nunawading Victoria 3131, Australia
E jobskin.production@smith-nephew.com

Smith & Nephew Pty Ltd Australia
85 Waterloo Road North Ryde NSW 2113 T +61 2 9857 3999

Smith & Nephew Ltd New Zealand
Unit A 36 Hillside Road Wairau Valley Aukland 0627 T +64 9 820 2840

Trademark of Smith & Nephew
SN12678 (05/2017)

Body Suit

Measurement Form – Page 1 of 3

Patient details (Please print
in black with CAPITAL letters)

Last name	
First name	Date



Sleeveless
Body Brief - **0530**



Body Brief
with Sleeves - **0531**



Singlet Body
Brief - **0556**



Singlet
Body Suit - **0557**
To distal measurement
above knee



Sleeveless
Body Suit - **0558**
To distal measurement
above knee




Body Suit
with Sleeves - **0560**
To distal measurement
above knee



Singlet Body Suit,
Full Leg - **0564B**



Sleeveless Body Suit,
Full Leg - **0562B**



Body Suit with Sleeves,
Full Leg - **0563B**

Options and modification			
Code	Qty	Description	
0034		Turtleneck - Powernet	<input type="checkbox"/> Attached <input type="checkbox"/> Detached <input type="text" value=""/> cm
0034		Turtleneck - Elastic	<input type="checkbox"/> Attached <input type="checkbox"/> Detached <input type="text" value=""/> 5cm <input type="text" value=""/> 2.5cm <input type="text" value=""/> 1.3cm
BRACUP		Bra Cup	Cup Size: <input type="text" value=""/>
1203		Postural Support (please indicate on TBSA Map)	
1203		Shoulder Bracing	<input type="checkbox"/> Left <input type="checkbox"/> Right
1190		Anterior Elbow Insert	<input type="checkbox"/> Left <input type="checkbox"/> Right
1175		Double Abdominal Panel	
1176		Uplift Panel	
1161		Soft Abdominal Panel (Pregnancy)	
1116		Ostomy Opening	
1177B		Knee Patch (Powernet)	<input type="checkbox"/> Left <input type="checkbox"/> Right
SBE		Silicone Backed Elastic	
0027		Padding - Soft (indicate on TBSA Map)	<input type="checkbox"/> Pocket <input type="checkbox"/> Attached
0027		Padding - High Density (indicate on TBSA Map)	<input type="checkbox"/> Pocket <input type="checkbox"/> Attached
1202		Motif	

Zippers			
Code	L	R	Description
1164			Arm Zippers (20 cm or less)
1165			Arm Zippers (over 20 cm)
1164			Leg Zippers (20 cm or less)
1165			Leg Zippers (over 20 cm)
ZCT			Zipper Velcro [®] Closure Tab <input type="checkbox"/> Arms <input type="checkbox"/> Legs <input type="checkbox"/> Body
1163B			Zipper Fastening Aids <input type="checkbox"/> Arms <input type="checkbox"/> Legs <input type="checkbox"/> Body
1172A			Body Zipper Loop (leather)
0027			Padded Zipper Lining <input type="checkbox"/> Arms <input type="checkbox"/> Legs <input type="checkbox"/> Body
0000			Satin Zipper Lining <input type="checkbox"/> Arms <input type="checkbox"/> Legs <input type="checkbox"/> Body
0000			Elastic Zipper Lining <input type="checkbox"/> Arms <input type="checkbox"/> Legs <input type="checkbox"/> Body

Linings			
Code	L	R	Description
0000			Lining – Axilla <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1167			Lining – Inside Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1168			Lining – Full Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1167B			Lining – Posterior Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
0040			Lining - Behind Knee <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1177			Lining - Anterior Knee <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1197			Fully Lined Garment <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1188			Lining - Elsewhere (indicate on TBSA map) <input type="checkbox"/> Satin <input type="checkbox"/> Cotton

Design choices			
Neck drop		Axillas	
<input type="checkbox"/> Front Neck Drop	<input type="text" value=""/> cm	Powernet	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> Back Neck Drop	<input type="text" value=""/> cm	Powernet + Lining	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> Scoop Front	<input type="checkbox"/> "V" Neck Front	Mesh	<input type="checkbox"/> Left <input type="checkbox"/> Right
<input type="checkbox"/> Wide Neck		Satin	<input type="checkbox"/> Left <input type="checkbox"/> Right
Breast Support		Closures	
<input type="checkbox"/> Bra Cups	<input type="checkbox"/> Bra Cup Lace	<input type="checkbox"/> Front Zipper	<input type="checkbox"/> No Velcro [®] Tab
<input type="checkbox"/> Contoured Breast		<input type="checkbox"/> Back Zipper	

Garment styles		
Code	Qty	Description
0530		Sleeveless Body Brief
0531		Body Brief with Sleeves
0556		Singlet Body Brief
0557		Singlet Body Suit
0558		Sleeveless Body Suit
0560		Body Suit with Sleeves
0564B		Singlet Body Suit, Full Leg
0562B		Sleeveless Body Suit, Full Leg
0563B		Body Suit with Sleeves, Full Leg

Pubis styles	
Male	Female
1. Open <input type="checkbox"/>	1. Powernet & Lining <input type="checkbox"/>
2. Diagonal Fly <input type="checkbox"/>	2. Satin <input type="checkbox"/>
3. Gusset (fewer inseams) <input type="checkbox"/>	3. Mesh <input type="checkbox"/>
4. Pouch <input type="checkbox"/>	4. Open <input type="checkbox"/>
Options for styles 2,3 & 4	
Powernet <input type="checkbox"/>	Velcro [®] <input type="checkbox"/>
Powernet & Lining <input type="checkbox"/>	Snaplock <input type="checkbox"/>
Satin <input type="checkbox"/>	

Body Suit

Measurement Form – Page 2 of 3


Patient details (Please print
in black with CAPITAL letters)

Last name

First name

Date

Measuring points



Measurements

cm / mm

Neck circumference	1				
Left shoulder circumference	2				
Right shoulder circumference	3				
Chest circumference	4				
Circumference over nipple line	5				
Circumference just under breast (level of breast attachment to chest)	6				
Waist circumference	7				
Circumference of largest part of belly	8				
Circumference of midpoint: Waist and Buttocks	9				
Circumference of largest part of buttocks	10				
Proximal Thigh Circumference: Left	11				
Proximal Thigh Circumference: Right	12				
C7 to nipple line	13				
C7 to just under breast	14				
C7 to waist	15				
C7 to largest part of belly	16				
C7 to largest part of buttocks	17				
C7 to fold of buttocks	18				
Side seam to side seam under axilla on back	19				
Side seam to side seam at waist on back	20				
Gluteal fold to end of support	21				
Shoulder width	22				
Singlet/sleeveless garment shoulder width	22a				
Top arm tape to neck (arm down by side)	23				
Under breast to waist on front	24				

Patient to wear well fitted bra during
measurement of Breast Support.

