

MEASUREMENT  
FORM

 **smith&nephew**  
**JOBSKIN<sup>®</sup>**  
Made to Measure  
Compression Garments



# Glove, Gauntlet

SN12677 (05/2017)

Glove, Gauntlet

Order/Reorder Form

All measurements taken with no tension and in centimetres only.

Design options



Contrast Colour Inserts



Slant inserts



Suede grip



Thenar eminence seam-free



Gauntlet



Palmar reinforcement



Silicone grip



Finger stall with wrist attachment

Gloves & gauntlets definition, hints and tips

**Gauntlet measurements required:** palm, thumb, wrist and any arm measurements needed. The hand outline must be completed.


The hand drawing must be perpendicular to the edge of the paper, middle finger perpendicular to the wrist, and palm down. The fingers and palm should be flat against the paper. The important points to mark exactly are the web spaces, wrist and finishing point fingers.

If fingers are contracted measure from the web space along the side of the fingers and follow the curve to the desired length of each finger. **A glove to wrist** may extend up to 3 tape.

**Finger Circumferences:** Take with a JOBSKIN finger tape and measure in skin contact - not tightly.

**Open Tips:** Measure only the length desired from the Finger Web Space (palmar aspect).

**Slant Inserts:** Enable a snug fit on the dorsal surface and give extra pressure in the finger webs.




Use the finger tape to measure finger circumferences of hand. Use the **teal** measuring tape for the arm.

**Web Spacers:** Are designed to give extra pressure in the web spaces and are worn over gloves. They can be ordered with or without a thumb piece.

**Zippers:** May be used to allow for inserts, e.g. silicone gel sheets.

Arm measuring tips

With the arm straight and horizontal, proceed measuring using a **Moss** JOBSKIN arm tape.




**For a Full Length Sleeve** begin by securing tape 9 directly around the elbow. The last tape should fall just below the wrist bone (distal to styloid process). Pleating at the spine occurs **ONLY** at the wrist and axilla if required.

**For a Forearm Sleeve** begin tape at wrist and pleat at proximal end if required.


Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram).



Reading the tapes

Always read along the straight edge and in direction of arrows as indicated below.



Email, fax or return forms to:

**Smith & Nephew Pty Ltd Australia**  
Unit 1, Level 1, 100 Station Street  
Nunawading  
Victoria 3131, Australia  
T 1800 772 577  
F 1800 773 577  
www.smith-nephew.com/australia  
E Jobskin.cs@smith-nephew.com

**Smith & Nephew Ltd New Zealand**  
621 Rosebank Road, Avondale  
PO Box 422, Shortland Street  
Auckland 1140, New Zealand  
T 0800 765 627  
F 0800 465 278  
www.smith-nephew.com/new-zealand

Garment returns

**Smith & Nephew Pty Ltd Australia**  
Unit 1, Level 1, 100 Station Street  
Nunawading Victoria 3131, Australia  
E jobskin.production@smith-nephew.com

**Smith & Nephew Pty Ltd Australia**  
85 Waterloo Road North Ryde NSW 2113 T +61 2 9857 3999

**Smith & Nephew Ltd New Zealand**  
Unit A 36 Hillside Road Wairau Valley Aukland 0627 T +64 9 820 2840

®Trademark of Smith & Nephew  
SN12677 (05/2017)

Last name

First name

Date

### How to complete hand outline page



Gauntlet – MCPs to Wrist **0505**



Glove to Wrist **0535**



Interdigital Web Spacer **0536**



Mitten (For amputation at MCP level) **0537**



Gauntlet to Axilla **0502**



Glove to Axilla **0532**



Gauntlet to Elbow **0516**



Glove to Elbow **0534**




Gauntlet to Axilla with Shoulder Attachment **0504**



Glove to Axilla with Shoulder Attachment **0533**

CLEARLY MARK END OF GARMENT ON FINGERS



OT - Open Tip  
CT - Closed Tip

CLEARLY MARK  
WEB SPACES

CLEARLY MARK  
WRIST LINE

### Powernet colour selection guide

**Fabric & Thread:** Black, Brown, Beige, Light Beige, Purple, Navy, Royal, Light Blue, Red, Pink, White, Camouflage

**Contrast Thread:** Orange, Hot Pink, Lime, Yellow

Fabric	Thread	Qty
1	1	
2	2	
3	3	
Must provide extra colour options		
1	1	
2	2	

### Hand/arm zipper locations

	Location		Full Length	Or Specify Length cm
	L	R		
Dorsal				
Ulnar				
Palmar				
Other:				

### Comments

### Garment styles

Code	L	R	Description
0505			Gauntlet - MCPs to wrist
0535			Glove to Wrist <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0536			Interdigital Web Spacer <input type="checkbox"/> Thumb <input type="checkbox"/> No Thumb
0537			Mitten to Wrist (amputee)
0502			Gauntlet to Axilla
0532			Glove to Axilla <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0516			Gauntlet to Elbow
0534			Glove to Elbow <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0504			Gauntlet with Shoulder Att.
0533			Glove with Shoulder Att. <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips
0025			Individual Finger Stall <input type="checkbox"/> Open Tips <input type="checkbox"/> Closed Tips <input type="checkbox"/> Wrist Band

### Options and modifications

Code	L	R	Description
0021			Reinforced Palm (Powernet)
1183			Palm Grip <input type="checkbox"/> Silicone <input type="checkbox"/> Suede
1188			Positional Reinforcement (indicate on hand outline)
1169			Slant Inserts
0027			Padding - Soft (indicate on TBSA Map). <input type="checkbox"/> Pocket <input type="checkbox"/> Attached
0027			Padding - High Density (indicate on TBSA Map). <input type="checkbox"/> Pocket <input type="checkbox"/> Attached
1202			Motif
SBE			Silicone Backed Elastic
1190			Anterior Elbow Insert .....degrees
0000			Elbow Contracture Seam .....degrees
GCCINSERT			Contrast Colour Inserts (specify 3 choices) 1..... 2..... 3.....
0021			Lining-Palm
1167			Lining-Inside Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1168			Lining-Full Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1167B			Lining-Posterior Elbow <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1188			Lining (specify location) <input type="checkbox"/> Satin <input type="checkbox"/> Cotton
1195			Full Lined Garment <input type="checkbox"/> Soft Satin <input type="checkbox"/> Cotton
1164			Zippers (20cm or less)
1165			Zipper (over 20cm)
1172A			Zipper Loop (leather)
ZCT			Zipper Closure Velcro® Tab (top of zipper)
1163B			Zipper Fastening Aids (along length of zipper)
0000			Zipper Lining <input type="checkbox"/> Satin <input type="checkbox"/> Elastic
0027			Padded Zipper Lining

Patient details (Please print in black with CAPITAL letters)

Last name

First name

Date

Desired finger lengths

Diagram showing desired finger lengths for the left and right hands. The diagram includes boxes for measurements at the thumb, index, middle, ring, and little fingers, as well as the wrist and elbow. Arrows indicate the direction of measurement. White boxes are lengths and grey boxes are circumferences.

Arm circumferences

Left (cm/mm)	Wrist End	Right (cm/mm)	Pleat
	-3		
	-1½		
	0		
	1½		
	3		
	4½		
	6		
	7½		
	Elbow 9		
	10½		
	12		
	13½		
	15		
	16½		
	18		
	19½		

Diagram showing arm circumferences for the left and right arms. The diagram includes boxes for measurements at the wrist, elbow, and axilla. Arrows indicate the direction of measurement. White boxes are lengths and grey boxes are circumferences.

White boxes are lengths and grey boxes are circumferences

3cm past wrist

6cm past wrist

For garments finishing above elbow measure length from wrist to end of garment

Desired finger lengths

Diagram showing desired finger lengths for the left and right hands. The diagram includes boxes for measurements at the thumb, index, middle, ring, and little fingers, as well as the wrist and elbow. Arrows indicate the direction of measurement. White boxes are lengths and grey boxes are circumferences.

White boxes are lengths and grey boxes are circumferences

3cm past wrist

6cm past wrist

For garments finishing above elbow measure length from wrist to end of garment

Shoulder attachment measurements

1. Distance last arm tape to inside edge of bra strap	2. Top of shoulder, circumferentially under opposite breast	3. Top of shoulder, circumferentially under opposite axilla	4. Top of shoulder, diagonal to waist
Bra Loop (cm)	Female Body Strap (cm)	Male/Child's Body Strap (cm)	Waist Attachment (cm)
1	2	3	4
2a, 3a, 4a: Distance last arm tape to base of neck	2a	3a	4a

