

MEASUREMENT
FORM

 **smith&nephew**
JOBSKIN[®]
Made to Measure
Compression Garments



Head, Neck, Face Mask


SN12681 (05/2017)

Head and Neck Face Mask

Order/Reorder Form

All measurements taken with no tension and in centimetres only.

Design Options



Measuring tips

Face mask & chin strap instructions

- Use regular metric tape measure for circumferences.
- Use two rulers to measure 4. Hold one ruler horizontally under the chin. Hold one ruler vertically and at right angles to the horizontal ruler and read off length. Ensure child is safe to wear neck collar.

- Step 1:** Select garment style
- Step 2:** Measure patient following the numbered measurements under the picture of your chosen garment style
- Step 3:** Complete body measurement grid
- Include patient photos with order if possible
 - When multiple head and neck garments are required, initially order one only to ensure accurate fit, then reorder
 - It is preferable for first garment to not have ear holes cut, so the Therapist can mark correct position

Email, fax or return forms to:

**Smith & Nephew Pty Ltd
Australia**
Unit 1, Level 1, 100 Station Street
Nunawading
Victoria 3131, Australia
T 1800 772 577
F 1800 773 577
www.smith-nephew.com/australia
E Jobskin.cs@smith-nephew.com

**Smith & Nephew Ltd
New Zealand**
621 Rosebank Road, Avondale
PO Box 422, Shortland Street
Auckland 1140, New Zealand
T 0800 765 627
F 0800 465 278
www.smith-nephew.com/new-zealand

Garment returns

Smith & Nephew Pty Ltd Australia
Unit 1, Level 1, 100 Station Street
Nunawading Victoria 3131, Australia
E jobskin.production@smith-nephew.com

Smith & Nephew Pty Ltd Australia
85 Waterloo Road North Ryde NSW 2113 T +61 2 9857 3999

Smith & Nephew Ltd New Zealand
Unit A 36 Hillside Road Wairau Valley Aukland 0627 T +64 9 820 2840

®Trademark of Smith & Nephew
SN12681 (05/2017)

Head, Neck, Face Mask

Last name	First name	Date
-----------	------------	------

Garment styles			Options & modifications			Powernet colour selection guide		
Code	Qty	Description	Code	Qty	Description	Fabric & Thread: Beige, Light Beige, Black, Brown, Navy, Royal, Light Blue, Purple, Red, Pink, White, Camouflage Contrast Thread: Orange, Hot Pink, Lime, Yellow		
0550		Chin Strap	1188		Lining - Anywhere (Specify location)	Fabric	Thread	Qty
0549		Modified Chin Strap (extends behind ear)	LIPSTRAP		Lip Strap used with Open Face mask	1	1	
0546		Open Head and Open Face Mask	0543		Nose Covering	2	2	
1158		Open Face Mask	1188		Nose contouring (double layer of p/net, used with nose covering only)	3	3	
0541		Open Head Mask	1188		Ear Reinforcement <input type="checkbox"/> L <input type="checkbox"/> R	Must provide extra colour options		
0540		Face Mask	1188		Eye Reinforcement <input type="checkbox"/> L <input type="checkbox"/> R	1	1	
0545		Chin Extension Collar	1188		Ear Pressure Flap <input type="checkbox"/> L <input type="checkbox"/> R	2	2	
0034		Attached Turtleneck - Elastic <input type="checkbox"/> 5cm <input type="checkbox"/> 2.5cm <input type="checkbox"/> 1.3cm	1202		Motif			
0034		Detached Turtleneck - Powernet Height <input type="checkbox"/> <input type="checkbox"/> cm	0027		Padding			
			Cut Earholes <input type="checkbox"/> Yes					
						Mask closure		

[illegible]

Comments	

