

MEASUREMENT FORM

MEASUREMENT
FORM

 **smith&nephew**
JOBSKIN[®]
Made to Measure
Compression Garments



Knee High, Anklet, Toe Piece

SN12675 (05/2017)

Knee High, Anklet, Toe Piece

Order/Reorder Form

All measurements taken with no tension and in centimetres only.

Design options




Foot outline measuring tips

Centre the patient's foot on the foot outline page as marked. Ensure that the patient is standing or sitting with the foot flat on the ground. Mark the following clearly:


- Back of the heel
- Webspaces
- Toe lengths (if measuring a toepiece or glove)
- Where the garment will finish

Use the grey measuring tape for these garments. Measure in centimetres only.




Pleat instructions

When pleating the last tape, measure the distance between the pleated tape and its neighbour from mid point to mid point (see diagram). Measurements must be recorded in centimetres only.



Reading the tapes

Always read along the straight edge and in direction of arrows as indicated below.



Email, fax or return forms to:

Smith & Nephew Pty Ltd Australia
Unit 1, Level 1, 100 Station Street
Nunawading
Victoria 3131, Australia
T 1800 772 577
F 1800 773 577
www.smith-nephew.com/australia
E Jobskin.cs@smith-nephew.com

Smith & Nephew Ltd New Zealand
621 Rosebank Road, Avondale
PO Box 422, Shortland Street
Auckland 1140, New Zealand
T 0800 765 627
F 0800 465 278
www.smith-nephew.com/new-zealand

Garment returns

Smith & Nephew Pty Ltd Australia
Unit 1, Level 1, 100 Station Street
Nunawading Victoria 3131, Australia
E jobskin.production@smith-nephew.com

Smith & Nephew Pty Ltd Australia
85 Waterloo Road North Ryde NSW 2113 T +61 2 9857 3999

Smith & Nephew Ltd New Zealand
Unit A 36 Hillside Road Wairau Valley Aukland 0627 T +64 9 820 2840

®Trademark of Smith & Nephew
SN12675 (05/2017)

Patient details (Please print
in black with CAPITAL letters)

Last name

First name

Date

Toe circumference
LEFT

☐ Open Tips
☐ Closed Tips

Toe circumference
RIGHT

☐ Open Tips
☐ Closed Tips

Diagram illustrating the measurement process for the Toe Piece, showing the placement of the piece on the foot and the corresponding measurement points (A and B) for both the left and right feet.

Labels for the diagram include:

- Desired toe lengths
- 1, 2, 3, 4, 5 (Measurement points)
- A, B (Measurement points)
- Length
- Circumference of MTP (for toe piece only)
- Circumference at end of toe piece

Instructions for measurement:

- Toe circumference 'B' should be 'base of toe' & 'A' should be 'at point of finished length' (mm)
- (White boxes are lengths and grey boxes are circumferences)

Note: FOOT OUTLINE (SEE: PG 3) MUST ACCOMPANY THIS FORM

Patient details (Please print
in black with CAPITAL letters)

Last name

First name

Date



JOBSKIN®

Made to Measure
Compression Garments

27	26	25	24	23	22	21	20	19	18	17	16	15	14	13	12	11	10	9	8	7	6	5	4	3	2	1	0	cm
For open tips - indicate on foot outline desired length																												
Please ensure all bases of webspaces are clearly marked																												
Place third toe on line																												
Heel End																												
Toe end																												